

Differential Micrometer Heads (Maintenance parts) Nut Drivers

MHF
MNS

MHF (Travel Coarse: 0 – 13mm, Fine: 0 – 0.2mm)





Catalog Code **W8003**



MHF2-13

MHF2-13F



Attention

- ▶ Depending on the stage type, there is the possibility of shortened travel or interference with mounting parts. Contact our Sales Division when you swapping micrometer heads.

Specifications

Part Number	Travel [mm]	Type	Readable Resolution [mm]	Lead of Actuator [mm]	Load Capacity (Static load) [N]	Weight [kg]
MHF2-13	Coarse: 0 – 13 Fine: 0 – 0.2	Coarse/Fine	Coarse: 0.01 Fine: 0.0005	Coarse: 0.5 Fine: 0.025	9.8 (1kgf)	0.1
MHF2-13F	Coarse: 0 – 13 Fine: 0 – 0.2	Coarse/Fine	Coarse: 0.01 Fine: 0.0005	Coarse: 0.5 Fine: 0.025	9.8 (1kgf)	0.1

MNS

Catalog Code **W8004**

Special tool required to remove or secure micrometer heads.

- Special nuts are used for fastening micrometers on some manual stages (TSD and TSDS series). This custom nut driver makes removing and attaching these nuts easier and faster.



Guide

- ▶ This product comes as a set of two head sizes. Refer to the following table for compatible stages.

Head Size [mm]	Part Number of Compatible Stage
For φ8	MDC-6.5/6.5R
For φ12.5	TSD-40*/60/65/-80, TSDC-40/60, TSDH-40/60
For φ14.5	TSD-100/TSD-120

* TSD-40 C type only.

Attention

- ▶ Not all manual stages (TSD and TSDS series) are compatible. Contact our Sales Division if you have any questions.
- ▶ Micrometer heads can be swapped with a finepitch screw upon request. Note self installation may not be covered by the warranty of the stage's warranty.


Specifications

Part Number	MNS-8
	MNS-12.5-14.5

Outline Drawing

MNS-12.5-14.5

- Head for φ12.5 (mounted on the grip) (compatible with TSD-40/60/65/-80, TSDC-40/60, TSDH-40/60)



- Head for φ14.5 (compatible with TSD-100 and 120)

