

April, 2021

Dear Valued Customers,

SIGMAKOKI CO., LTD.

Sales Department

**Notice of change of the specs and the part number of
2 axis stage controller (SHOT-702)**

Thank you very much for your constant patronage.

Due to the end of production of the main internal parts of the **2-axis stage controller (SHOT-702)**, some specifications and part numbers will be changed.

SHOT-702 is scheduled to be discontinued in the fall of 2021.

Sales of the successor model, **SHOT-702H**, are scheduled to commence around May 2021.

Please refer to the table below for the differences in specifications.

Should you have any questions or concerns, please do not hesitate to contact us.
Thank you for your continued support.

2 axis stage controller Catalogue code : W9045

| | End of sales | Successor product |
|---------------------------------------|---------------------|--|
| Part Number | SHOT-702 | SHOT-702H |
| Supply Voltage | AC100-240V | DC24V 3.4A※ |
| Power consumption | 50VA | 81.6VA |
| Drive current (RUN) (A/phase) | 0.1~1.1 | 0.25~1.4 |
| Stop current (STOP)(A/phase) | 0.1~0.7 | 50% of RUN (fixed) |
| Interface | RS232C | RS232C/ USB (virtual COM port) |
| Electromagnetic brake power supply | Not equipped | Equipped (STAGE connector 9,10 pins) |
| Dimensions(WxDxH) [mm] | 260x280x70 | 180x119x55 |
| Weight(kg) | 2.8 | 1.0 |

※Dedicated AC adapter AC100-240V included